

# Organic LED imaging by AP/MALDI(ng) UHR / Orbitrap HRMS

Dana El Assad and Gilles Frache

Materials Research and Technology Department, Luxembourg Institute of Science and Technology (LIST),  
 a MassTech European Application Lab

## INTRODUCTION

MassTech AP/MALDI(ng) UHR module is a flexible add-on ion source for existing Thermo LC/HRMS equipment. This Application Note presents AP-MALDI MS analyses and imaging of an OLED display based on a Thermo *LTQ/Orbitrap Elite* instrument.

## WORKFLOW

### Sample Preparation:

- 1) A broken smartphone was disassembled to access the OLED panel
- 2) A portion of the Display was cut and adapted on MassTech Glass Slide holder and adjusted with spacers to meet the standard working distance
- 3) The homogeneous protective layer was removed using sticky tape



### Data Acquisition:

- 1) Definition of image parameters (CSR mode, image dimension 500 x 500microns, 10 micron pixel size) in MassTech *Target* (a control software)  
 → Parameters are saved as .xml file
- 2) Molecular imaging using AP/MALDI(ng) UHR ion source (MassTech) with *LTQ/Orbitrap Elite* (Thermo) high resolution MS
- 3) Data handling with *ImageQuest* (Thermo) using native files (XML and RAW files)



Fig.1: MassTech AP/MALDI(ng) UHR add-on source

## KEY CHARACTERISTICS AND BENEFITS

- MALDI imaging capabilities down to 10 micron lateral resolution
- Switch from the AP-MALDI configuration to LCMS configuration within 5 minutes.
- Provides MALDI analysis and imaging capabilities to high-end LCMS instruments
- Laser focus can be adjusted from 10 to 30 microns

## Available modes of operation

- Sequential MALDI MS analysis of multiple sample spots on AB Sciex *OptiTOF 192* or Bruker *MTP96* MALDI target plates using spiral or raster motion with tunable parameters
- MALDI MS imaging using Pixel-Map mode, with tunable in-pixel motion for various tissue MS analyses (HRMS, SRM, etc).
- MALDI MS imaging using Constant Speed Raster (CSR) mode for increased pixel data acquisition rates

### Case Study:

AP-MALDI UHR imaging capability has been applied for OLED ingredients MS imaging using a Thermo HRMS equipment. High lateral resolution image was demonstrated using spot size below 10 micrometers to produce MALDI MS images over a field of view of 500  $\mu\text{m}$  x 500 $\mu\text{m}$  (Fig. 2).

### Ingredient identifications:

AP-MALDI UHR imaging allows for accurate mass measurement of ions produced from a 10 $\mu\text{m}$ x10 $\mu\text{m}$  spot size (within +/-3mDa without lockmass). Structural confirmation may be obtained from simultaneous in-pixel MS/MS analysis.




Fig.2: MassTech AP/MALDI(ng) UHR imaging of a portion of an OLED display (R:  $m/z=638.2685$ , G:  $m/z=729.2635$ , B:  $m/z=675.2890$ )




Fig.3: Differential Single scan High-Resolution mass spectra extracted from areas in Red (top), green (middle) and Blue (bottom) in Fig.2. based on background subtraction feature in XCalibur




Fig.4: Tentative assignment of  $m/z=638.270$  detected in areas in Red in Fig.2. (calculated MS for  $C_{48}H_{32}N_2^+$  (top) vs. experimental spectrum and its AP-MALDI HRMS image (bottom))

AP-MALDI HRMS Imaging reveals molecular contrasts from single 10 $\mu\text{m}$ x10 $\mu\text{m}$  areas. Application to organic microelectronics may be useful to study aging effects or failures.

MassTech offers a range of analytical ionization sources, accessories and complete instruments for advanced analytical platforms. An AP-MALDI UHR ion source is available worldwide exclusively from MassTech and authorized resellers, such as KR Analytical in Europe ([www.kranalytical.co.uk](http://www.kranalytical.co.uk)).

To request further information, please contact: [sales@apmaldi.com](mailto:sales@apmaldi.com)

MassTech selected LIST as European Application Lab: [www.list.lu](http://www.list.lu)